



Plate Carriers

Stainless Steel Carriers for Easy Transportation of Plates

Environmental monitoring in cleanrooms is a critical task used to demonstrate compliance to standards of cleanliness. It is equally critical to maintain sample integrity until the final reporting of results several days later.

Developed and designed by Cherwell, the range of stainless steel Plate Carriers offer an easy way to handle and transport plates to and from sampling locations and the laboratory. Available for 90mm Petri dishes and Contact plates, each carrier can hold up to 20 plates.

Applications

The Plate Carriers can be used within the laboratory and when deploying environmental monitoring, the key uses are:

- **Routine clean room monitoring**
- **Plate storage within the laboratory**
- **Segregation during incubation**

Features & Benefits

Stainless steel construction

Each carrier is made completely from stainless steel to provide a cleanroom compatible finish that is easily cleaned and/or autoclaved.

Easy clean design

The carriers have no fixings or tight corners, making cleaning easy and effective.

Swing handle

Each carrier has a pair of swing handles which make the carriers easy to handle. The handles swing out of the way allowing unobstructed loading of the carrier.



Technical Data

Product Specifications - 55mm Contact Plate

Materials	304 Stainless Steel
Overall Dimensions (mm): *with handle folded down	145 x 85 x 135* (L x W x H)
Weight (kg):	0.5

Product Specifications - 90mm Petri Dish

Materials	304 Stainless Steel
Overall Dimensions (mm): *with handle folded down	194 x 107 x 160* (L x W x H)
Weight (kg):	0.65

Ordering Information

Part Number	Description	Packaging
3.2003	Plate Carrier S/S – Contact (20)	Each
3.2004	Plate Carrier S/S – Petri (20)	Each

Part Number	Description	Packaging
3.2001	Settle Plate Stand S/S - single	Each
3.2006	Settle Plate Stand S/S - twin	Each

Terms & Conditions

All Plate Carriers are supplied with a 12 months warranty. We reserve the right to alter specifications without prior notice. E&OE

cherwell 
an analytichem company

7 & 8 Launton Business Centre, Murdock Road, Bicester, OX26 4XB, United Kingdom
T. +44 (0)1869 355500 F. +44 (0)1869 355545 E. sales@cherwell-labs.co.uk

www.cherwell-labs.co.uk